



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

EMULSIFIER / HOMOGENIZER CE



OT-E-54

Fitted with motor having speed regulator with maximum speed upto 4000 RPM with perforated S.S. Disc (Complete on stand).

Model No. Available with motor size:

OT-E-54/A	(A) 1/20 H.P. Motor
OT-E-54/B	(B) 1/4 H.P. Heavy Duty Motor
OT-E-54/C	(C) 1/2 H.P. Heavy Duty Motor
OT-E-54/D	(D) 2/3 H.P. Heavy Duty Motor
OT-E-54/E	(E) 1 H.P. Heavy Duty Motor

Optional :

- 1) Digital RPM Indicator
- 2) Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@optics technology.in,

optics technology@gmail.com,

Website : www.optics technology.in