



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

GEL ELECTROPHORESIS UNIT

HORIZONTAL / VERTICAL

CE



OT-GEU-231

Electrophoresis of nucleic acid (DNA/RNA) in agarose gel has become an essential step in every walk of biotechnology today. High buffer-gel ratio is kept in all our systems whereby current flow is minimized; temperature rise and pH, difference are also kept at minimum levels.

This results in fast and good resolution of bands.

FEATURES:

- Basic unit with platinum electrode assembly.
- Comb 2mm thick
- 10 well and thick comb stand
- Connecting cord

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in