



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

GYRATORY SIEVE SHAKER Heavy Duty CE



OT-GSS-112

**GYRATORY
SIEVE SHAKER**

(WITH MOTOR & MECHANICAL TIMER)

(For Size 8" Sieve)

Produces upward and downward movement but is inclined from the vertical axis and the direction of inclination changes progressively in clockwise direction.

Optional:

- Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in