



# OPTICS TECHNOLOGY

SINCE 1978

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## LAB CENTRIFUGE

CE



OT-LC-245

### SPECIAL FEATURES :

Easily operated electronic lid lock, Speed regulator, safety Fuse.  
Elegant front panel.

With a wide choice of rotor heads and adaptors this unit is truly versatile. Digital RPM Indicator and Digital Timer.

Supplied with Rotar Head having 8 x 15 ml

### TECHNICAL DATA :

Max. Speed : 4000 RPM

OR

5000 RPM

**Note:** Centrifuges with different speed ranges and cooling models  
and with different rotors also available on request.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)