



LEAK TEST APPARATUS CE

(GMP MODEL)



OT-LTA-25

It's the most ideal 'LEAK TEST' instrument for checking Blister Packs and Strips / Pouches. The Vacuum Pump ensures that a high level of vacuum is created in the shortest possible time. A vacuum gauge is provided to indicate the vacuum level. A 0 to 5 minute timer helps the operator attend to other chores at a time. It makes the conventional desiccator with a huge vacuum pump, obsolete. And the die - pressed. Polycarbonate vacuum desiccator sustains the vacuum for a long time. This is connected to the isolation valve which helps disconnect the source, thereby avoiding the need for running the vacuum pump continuously.

Specifications :	OT-LTA-25/10	OT-LTA-25/12
Desiccator Diameter :	250mm (10" Dia)	300mm (12" Dia)
Rating :	5A 230V Single Phase	
Vacuum Level :	30" Hg	
Timer :	15 minute (Mechanical/Digital)	
Desiccator :	Polycarbonate	
Body :	Stainless Steel 304	
Actual Size (Inches) :	20" (W) x 14" (D) x 10.5" (H)	

Optional:

1. Digital Timer.
2. Microprocessor based Touch screen Display system with Pen Drive for data logging, printing through PC Communication & saving of data in PC