



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MICRO TISSUE HOMOGENIZER CE



This lightweight micro stirrer is ideal for use with tissue grinders having stainless steel shaft. Supplied complete with a stirrer/controller, tripod with a vertical support rod, clamp holder and a coupling Teflon paddle.

With Tissue Flask (Max. Speed 4,000-10,000 RPM with variable speed regulator).

Model No. Capacity

OT-MTH-172-A 10 ml

OT-MTH-172-B 30 ml

OT-MTH-172-C 50 ml

OT-MTH-172-D 100 ml

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in