



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MICROBIOLOGICAL BLENDER CE

(WITH AUTOCLAVABLE S.S. JAR)



OT-MBB-30

Technical Details:

- RPM** : 14000 - 18000 RPM Approx.
- Watts** : 450 Watts
- Timer** : Selective Timer Cycle from 1 to 5 minutes for automatic blending or manual control setting for continuous running.
- Base** : Aluminium Cast - Baked Enamel Finish
- Container** : Stainless Steel Jar with S.S. Blades & S.S. Cap.
- Capacity** : (a) 250ml (b) 500ml (c) 1000ml.

Model No. :

OT-MBB-30/250 : 250ml containers of Stainless Steel

OT-MBB-30/500 : 500ml containers of Stainless Steel

OT-MBB-30/1000 : 1000ml containers of Stainless Steel

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in