



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



QUARTZ DOUBLE DISTILLATION CE



OT-QDS-125

Quartz Double Distiller - with Replaceable Silica Heater and ripple Built-in Safety.

“Optics Technology” unit consist of main boiler with a inbuilt heater enclosed in quartz glass tube with auto water lever adjuster and bottom outlet for cleaning waste discharged the whole unit is mounted on epoxy coated metal box stand with inbuilt power supply unit.

Specifications:

- Replaceable silica Heater which eliminates frequent shutdown of the equipment.
- Gate valve which conserves cooling water wastage during Power failure.
- Flow switch which shuts off the unit if cooling water to condenser fails.
- Pressure temperature switch (PTS) which switch off the unit if the water in the boiler falls below the heater.
- Quartz Teflon stopcock for easy draining and cleaning of boiler.
- Quartz threaded connection for easy assembly and maintenance.
- Borosilicate Condenser & Quartz Boiler (Specific Conductivity: 08 - 1.0 $\mu\text{S}/\text{cm}$).

Capacity:

1.5 Liters

3 Liters

Optional :

Safety Control (Cut off Device)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in