



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

QUARTZ SINGLE DISTILLATION CE



OT-QSD-124

The unit is designed to produce a distillate of consistent purity the units continue to be most effective and economical method of producing ultra pure water for laboratory use.

Specification :

Voltage	: 220v
Biological Activity	: Pyrogen free
Heater	: Quartz
Boiler	: Quartz
Condenser	: Borosilicate

Capacity:

1.5 Litres
3 Litres

Optional :

Safety Control (Cut off Device)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in