



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

SOXHLET APPARATUS

CE



OT-SOX-43 MS

M.S. BODY (3 TEST)



OT-SOX-43 SS

S.S. BODY (6 TEST)

The complete soxhlet apparatus consists of solvent reservoir extraction assembly with 3/6 extraction places of various diameter to suit 100ml flasks to 250ml flasks and built in steam generator. Complete with stand, clamp holder soxhlet extraction (condenser) glass and extraction flask of 60 ml volume. Temp. range 50°C to 350°C. controlled by Sunvic energy regulator. Body is made of Mild Steel Duly Powder Coated/S.S. 304.

Model No.	Available in following models:
OT-SOX-43/3-MS/SS	3 Test Mantle Type OR 3 Test Hot Plate Type
OT-SOX-43/6-MS/SS	6 Test Mantle Type OR 6 Test Hot Plate Type

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in