



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

TABLET HARDNESS TESTER (PFIZER TYPE)



OT-THT-177 P

PFIZER TYPE

A handy, tough and direct dial reading type tablet hardness tester, working as pair of pliers, generates breaking force in single grip. The instrument offers accurate and reproducible result much faster other conventional hardness tester.

Specification :

Accuracy	: ± 0.2 kg
Least Count	: 0.2 kg
Safe Testing Force	: 20 kg
Maximum Tablet Diameter	: 25 mm

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in