



ULTRASONIC BATH (SONICATOR) CE



OT-UBS-168



OT-UBS-168

(S.S. 304 Body) With Soft touch Switches & Digital Timer

The offered system is widely used for cleaning different products like jewellery, lenses, watches, optical parts, surgical instruments and number of electronic equipments. Our offered system cleans all the products with the help of ultrasound using appropriate cleaning solvent.

Sonicator Bath Ultrasonic cleaners are especially effective for removing chemical residues from complex parts with crevices, small openings and joints. In addition to cleaning validation, ultrasonic cleaners and baths are used by chemical and pharmaceutical laboratories for solvent degassing and sample preparation.

Features:

- Noiseless operation
- Easy to use
- Low maintenance
- Ultrasonic frequency 40 KHZ \pm 3

Model No	Capacity	Tank size (L x B x H)	Wattage
OT-UBS-168-A	1.5 Ltr.	240 x 135 x 65 mm	60 W
OT-UBS-168-B	3.5 Ltr.	300 x 150 x 100 mm	120 W
OT-UBS-168-C	6.5 Ltr.	300 x 150 x 150 mm	200 W
OT-UBS-168-D	9 Ltr.	300 x 200 x 150 mm	250 W
OT-UBS-168-E	10 Ltr.	300 x 225 x 150 mm	250 W
OT-UBS-168-F	14 Ltr.	525 x 175 x 150 mm	300 W
OT-UBS-168-G	20 Ltr.	350 x 300 x 200 mm	400 W
OT-UBS-168-H	25 Ltr.	400 x 400 x 150 mm	500 W
OT-UBS-168-I	30 Ltr.	375 x 375 x 200 mm	600 W
OT-UBS-168-J	50 Ltr.	550 x 300 x 300 mm	750 W

Optional :

Digital Temperature Controller