



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

U.V. CABINET

CE



OT-UVC-34 MS

MS BODY



OT-UVC-34 SS

S.S 304 BODY



OT-UVC-34

Acrylic body

It is a portable darkroom design with a range of 3200-4000 A.U. for viewing paper & inspection chromatographs. It has a Wavelength of 254nm (Short wave) & 365nm (Longwave) and is suitable for other gross studies of fluorescence effects.

Chamber Size: 330 x 220 x 235 mm

Models:

OT-UVC-34-MS/SS-A	U.V. CABINET with Electronic Arrangement with Long, Short Wave (2 Tube)
OT-UVC-34-MS/SS-B	U.V. CABINET with Electronic Arrangement with Long & Short Wave and Daylight (3 Tube)
OT-UVC-34-C	U.V. CABINET with acrylic sheet shock Proof Body Short & Long Wave (2 Tube)
OT-UVC-34-D	U.V. CABINET with acrylic sheet shock Proof Body Short & Long Wave and Daylight (3 Tube)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in