



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

WATER BATH (SEROLOGICAL) CE

THERMOSTATICALLY CONTROLLED



OT-WBSE-37-MS

THERMOSTATICALLY
CONTROL



OT-WBSE-37-SS

INNER & OUTER S.S. 304
WITH THERMOSTAT

(Double walled) Inner Wall made up of S.S. 304 & Outer body made of M.S. Duly powder coated sheet/S.S. 304. The bath consists of 2 pilot lamps & on-off switch, immersion heaters. Temperature controlled by Thermostat from ambient temp. +5°C to 95°C $\pm$ 1°C OR Temp. can be controlled by a Microprocessor based PID Controller, ambient +5°C to 95°C $\pm$ 0.5°C.

Model No.	Size	Ltrs. (Approx.)
OT-WBSE-37-14 MS/SS	300 x 250 x 175 mm	14
OT-WBSE-37-18 MS/SS	350 x 250 x 200 mm	18
OT-WBSE-37-19 MS/SS	300 x 300 x 200 mm	19
OT-WBSE-37-22 MS/SS	350 x 300 x 200 mm	22
OT-WBSE-37-28 MS/SS	300 x 300 x 300 mm	28
OT-WBSE-37-28 MS/SS	450 x 300 x 200 mm	28
OT-WBSE-37-38 MS/SS	600 x 300 x 200 mm	38
OT-WBSE-37-48 MS/SS	600 x 400 x 200 mm	50
OT-WBSE-37-70 MS/SS	600 x 400 x 300 mm	75

- Optional :**
1. Microprocessor based P.I.D. Controller
 2. Water level cutoff device
 3. Circulation System/Stirrer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in