



# OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## ABBE REFRACTOMETER



Sugar 0 to 95% with both refractive index Sugar Scale.

It is a standard instrument for measuring the refractive index and dispersion of liquids and solids and determining the percentage contents of sugar in water solutions. Prism Box is water jacketed for circulating water at a constant temperature, if required, for precise measurement.

Measuring range	: nD 1.3000 to 1.7000
Measurement Accuracy	: ±0.1%

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com), Website : [www.opticstechnology.in](http://www.opticstechnology.in)