



BACTERIOLOGICAL INCUBATOR (LAB INCUBATOR)

CE



OT-BI-56 MS

INNER S.S. &
OUTER M.S. BODY



OT-BI-56 SS

INNER &
OUTER S.S. 304

Double walled, Inner chamber made of Stainless Steel 304. Outside Mild Steel painted in epoxy powder coating / stainless steel 304 with wire mesh adjustable shelves, two indicators and on /off switch to work on 220/230 volts A.C. Heating elements are placed in ribs at the bottom and sides. Temperature range Ambient Temp. **+5°C to 70°C**. The door has synthetic rubber gasket. With one number Air circulating Fan & Microprocessor based PID temperature controller. Specially Designed to maintain 37°C temp.

Model No.	Size Inside Chamber W x D x H	Ltrs. (App.)	No. of Trays
OT-BI-56/28-MS/SS	300 x 300 x 300 mm (12"x12"x12")	28	1
OT-BI-56/45-MS/SS	355 x 355 x 355 mm (14"x14"x14")	45	1
OT-BI-56/96-MS/SS	450 x 450 x 450 mm (18"x18"x18")	96	2
OT-BI-56/128-MS/SS	450 x 450 x 600 mm (18"x18"x24")	128	2
OT-BI-56/227-MS/SS	600 x 600 x 600 mm (24"x24"x24")	227	2
OT-BI-56/255-MS/SS	450 x 600 x 900 mm (24"x18"x36")	255	3
OT-BI-56/340-MS/SS	600 x 600 x 900 mm (24"x24"x36")	340	3
OT-BI-56/453-MS/SS	600 x 600 x 1200 mm (24"x24"x48")	453	4

OR ANY NON STANDARD SIZE

Optional:

- (A) Special heavy duty air circulating fan system.
- (B) Microprocessor based touch screen display system with Pen Drive for data logging, printing through PC communication & saving of data in PC.
- (C) Online transfer of data without any cable/wire Interference
- (D) 21 CFR software (for non-editable data)



OT-BI-56 MS