



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

B.O.D. INCUBATOR (DOUBLE DECKER)



OT-BOD-59 DD-MS

**DOUBLE DECKER
B.O.D. INCUBATOR**

(One above the other with individual temperature controlling temp. range + 5°C to 50°C)

Inner chamber made of S.S. 304 and outer body made of M.S. Duly powder coated body or stainless steel (S.S. 304) with independent controlling of temperature having independent control boxes and independent doors. Working size as below in each chamber both 2 Nos. compressors are fitted at the bottom of whole unit castor wheel mounted cabinet for easy mobility with 2 nos. microprocessor based PID controller.

Model No.	Inner Chamber Size W x D x H			Ltr.
OT-BOD-59DD-MS/SS-A	4 CU. Ft.	455 x 410 x 610mm		113 Ltr.
OT-BOD-59DD-MS/SS-B	6 CU. Ft.	505 x 415 x 830mm		173 Ltr.
OT-BOD-59DD-MS/SS-C	10 CU. Ft.	570 x 550 x 875mm		290 Ltr.

OR ANY SPECIAL SIZE

Optional:

1. 2 Nos. Microprocessor based Touch Screen Display System
2. Online transfer of data without any Cable/Wire interference
3. 21 CFR Software (For Non Editable Data)



2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in