



SINCE 1978

# OPTICS TECHNOLOGY

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



## BULK DENSITY APPARATUS

CE



OT-BDA-179

- This instrument is used to check the density of powder.
- Can Hold 2 No. cylinders, capacity : 100 ml.
- Based on and tapping mechanism, User has to select the No. of stokes in the counter
- Three models are available as given below Meet USP Standards.

### Model No.

**OT-BDA-179-T** Bulk Density Apparatus With **Timer**

**OT-BDA-179-D** Bulk Density Apparatus with **Digital Display of strokes**

**OT-BDA-179-A** **Automatic** bulk density apparatus

**BODY IS MADE OF MILD STEEL DULY POWDER COATED / S.S. 304**

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)