



# OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## BUTYRO REFRACTOMETER

(SUGAR AND OIL REFRACTOMETER)



OT-BR-14

This instrument is used to check the B.R. VALUE in Oil.

Divisions on the oil scale are sub-divided by Micrometer screw error due to temperature changes in sugar solutions can be eliminated by temperature corrections device provided. It is an indispensable tool for all industries using sugar, oil & many other similar substances. Butyro Refractometer Index has a Measuring range: 1.42 to 1.49 to an accuracy of 0.0002. Micrometer allows the accuracy of 0.01% on the oil (Butyro) scale.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in),

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com),

Website : [www.opticstechnology.in](http://www.opticstechnology.in)