



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL COLONY COUNTER

CE



OT-CC-31

REGULAR MODEL

The improved colony counter reveals colonies of bacteria clearly against a dark background. The illumination is brighter and completely uniform.

This principle increases the contrast between the colonies and the culture medium. It improves the visibility so those tiny colonies can be seen. As a result, the observer works efficiently and obtains accurate data quickly. Body is made of M.S./S.S. 304.

OT-CC-31-3-MS/SS Digital Colony Counter (3 Digits)

OT-CC-31-4-MS/SS Digital Colony Counter (4 Digits)

OT-CC-31-6-MS/SS Digital Colony Counter (6 Digits)

S.S. 304 BODY AS PER GMP IS ALSO AVAILABLE

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@optics technology.in,

optics technology@gmail.com,

Website : www.optics technology.in