



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL MELTING POINT APPARATUS

CE



OT-MPA-151-D

Shows digital display of Temperature in °C.

- For Determination of Melting Points upto 350°C
- Provision for 3 melting point tubes
- Magnifying lens and shadow less illumination of capillaries
- Fine temp. control
- It Has the facility to hold the temp.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in