



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DYNAMIC PASS BOX CE



OT-DPB-103 SS

HEPA Filter	: 0.3 μ down to 99.97%
Pre Filter	: 10 μ
Statics Pressure Manometer	
U.V. Light.	
Statically Balanced motor blower Assembly (Heavy Duty)	
Electronic Interlock System	
Very Less Electrical Load	
Energy Efficient Blowers	

Body is made of

Wooden / M.S. Duly Powder Coated / Stainless Steel (S.S. 304)

Model No.

OUTER SIZE: W x D x H

OT-DPB-103-W/MS/SS-A	18" x 18" x 18" (450 x 450 x 450mm)
OT-DPB-103-W/MS/SS-B	24" x 24" x 24" (600 x 600 x 600mm)

- Optional :** 1. Magnehelic Gauge.
2. LCD Display with soft touch switches

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in