



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DYNAMIC PASS BOX CE



OT-DPB-103 SS

HEPA Filter : 0.3 μ down to 99.97%

Pre Filter : 10 μ

Statics Pressure Manometer

U.V. Light.

Statically Balanced Motor Blower Assembly (Heavy Duty)

Electronic Interlock System

Very Less Electrical Load

Energy Efficient Blowers

Body is made of
Wooden / M.S. Duly Powder Coated / Stainless Steel (S.S. 304)

Model No.

OUTER SIZE:
W x D x H

OT-DPB-103- W/MS/SS-A 18" x 18" x 18" (450 x 450 x 450mm)

OT-DPB-103- W/MS/SS-B 24" x 24" x 24" (600 x 600 x 600mm)

Optional : 1. Magnehelic Gauge.
2. LCD Display with soft touch switches

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in