



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

FUME HOOD (SINGLE WALLED) CE



OT-FHSW-92 MS



BASE STORAGE
CABINET

OT-FHSW-92 SS

FUME HOOD (S.S.)



BASE STORAGE
CABINET

OT-FHSW-92 MS

M.S. BODY



BASE STORAGE
CABINET

OT-FHSW-92 SS

S.S. 304 BODY

Laboratory fume hoods are one of the most important components used to protect laboratory personnel from exposure to hazardous chemicals and agents used in the laboratory. Cabinet is made of :

- Wooden Body / M.S. Duly Powder Coated / Stainless Steel (304).
- Interiors acid resistant material coated with fire proof epoxy paint & further fixed on storage cabinet.
- Working table is covered with Granite
- Front door is made of 6mm clear Toughened glass.
- Counter weight balance Front Door.
- Vertical sliding facility & the front door can be stopped at any point.
- Exhaust blower consists of single phase motor connected with a PVC Exhaust Ducting. Standard ducting upto 6' is provided.

Model No. Working Size

OT-FHSW-92 W/MS/SS-A 2' x 2' x 2'

OT-FHSW-92 W/MS/SS-B 3' x 2' x 2'

OT-FHSW-92 W/MS/SS-C 4' x 2' x 2'

OT-FHSW-92 W/MS/SS-D 6' x 2' x 2'

OT-FHSW-92 W/MS/SS-E 8' x 2' x 2'

OR ANY NON STANDARD SIZE.

Optional:

1. Base Storage Cabinet.
2. P.P. Sink with Swan Neck Water Tap
3. Air flow monitor (LCD Display)
4. Automatic Motorized front door lifting arrangement
5. Touch Screen Display with Air Flow Monitor
6. Scrubber System



2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in