



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

HEATING MANTLE

CE

(With Energy Regulator)



OT-HM-40 MS

(M.S. BODY)



OT-HM-40 SS

(S.S. BODY)

Most suitable for labs, pharmaceuticals, chemical glass plants etc. Useful for heating inflammable liquids electric heating net is hand knitted from glass yarn (without any joint) and for durability & stands. Temperature up to 350°C. The body is made of Mild Steel duly powder coated / S.S. 304. Mantles lagged with special grade mineral wool. Most suitable to work on 220 volts A.C. Supply.

Model No.	Capacity
OT-HM-40 MS/SS 250	250ml
OT-HM-40 MS/SS 500	500ml
OT-HM-40 MS/SS 1000	1000ml
OT-HM-40 MS/SS 2000	2000ml
OT-HM-40 MS/SS 3000	3000ml
OT-HM-40 MS/SS 5000	5000ml
OT-HM-40 MS/SS 10000	10000ml
OT-HM-40 MS/SS 20000	20000ml

MOC: MS-Mild steel duly powder coated SS stainless steel kindly select the one required in model number.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in