



HOT AIR OVEN CE

(WITH TOUCH SCREEN DISPLAY & DATA LOGGING)



OT-HAO-55 MS-TS



OT-HAO-55 SS-TS

Double walled, Inner chamber made of Stainless Steel 304. Outside Mild Steel painted in epoxy powder coating / stainless steel 304 with wiremesh adjustable shelves, two indicators and on / off Switch. Heating elements are placed in ribs at the bottom and sides. Temperature range Ambient **+5°C to 250°C**. To work on 220 / 230 volts A.C. With one number Air circulating Fan.

- Touch Screen Display Menu driven settings for Controlling parameters Setting / RTC accurate Date Time Clock.
- Pen Drive for saving of data for printing through PC communication.

Model No.	Size	Ltrs. (App.)	No. of Trays
	W x D x H		
OT-HAO-55/28-MS/SS-TS	300 x 300 x 300 mm (12" x 12" x 12")	28	1
OT-HAO-55/45-MS/SS-TS	355 x 355 x 355 mm (14" x 14" x 14")	45	1
OT-HAO-55/96-MS/SS-TS	450 x 450 x 450 mm (18" x 18" x 18")	96	2
OT-HAO-55/128-MS/SS-TS	450 x 450 x 600 mm (18" x 18" x 24")	128	2
OT-HAO-55/227-MS/SS-TS	600 x 600 x 600 mm (24" x 24" x 24")	227	3
OT-HAO-55/255-MS/SS-TS	600 x 450 x 900 mm (24" x 18" x 36")	255	3
OT-HAO-55/340-MS/SS-TS	600 x 600 x 900 mm (24" x 24" x 36")	340	3
OT-HAO-55/453-MS/SS-TS	600 x 600 x 1200 mm (24" x 24" x 48")	453	4

Optional:

1. Special heavy duty air circulating fan system.
2. Online transfer of data without any cable/wire Interference
3. 21 CFR software (for non-editable data)

Calibration
is also Available

Note : Hot Air Oven is also available in high temperature (300°C)
(Model No. OT-HTO-55)