



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED
An ISO 10002:2018 CERTIFIED
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010
BS EN 12469:2000
ISO 29463-5:2022

HOT PLATE

(RECTANGULAR)



M.S. Body Cast Iron Top
with Energy Regulator

OT-HP-39 G (ER)



M.S. Body with S.S. Top
with Energy Regulator

OT-HP-39 H (ER)



S.S. Body with S.S. Top
with Energy Regulator

OT-HP-39 I (ER)

Technical Details:

- Max. Temp. upto 350°C
- With Energy Regulator
- Plate has insulation
- Elements To work on 220/230 Volts A.C.

MOC (Material of Construction):

- (a) M.S. Body with Cast Iron Top
- (b) M.S. Body with S.S. Top 316
- (c) S.S. Body with S.S. Top 316

Working Size:

- (A) 10" x 12" x 6" (250 x 300 x 150 mm)
- (B) 10" x 16" x 6" (250 x 400 x 150 mm)
- (C) 12" x 18" x 6" (300 x 450 x 150 mm)
- (D) 18" x 24" x 6" (450 x 600 x 150 mm)
- (E) 24" x 36" x 6" (600 x 900 x 150 mm)

Optional:

- Microprocessor Based PID Controller



M.S. Body Cast Iron Top
with PID Controller

OT-HP-39 J (PID)



M.S. Body with S.S. Top
with Energy Regulator

OT-HP-39 K (PID)



S.S. Body with S.S. Top
with PID Controller

OT-HP-39 L (PID)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in