



# OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## INCUBATOR SHAKER (Orbital Shaker Heating Only) CE



OT-INS-61 MS



OT-INS-61 MS

- With Variable speed and Microprocessor Based PID Controller
- Inner Stainless Steel & Outer M.S. / S.S. 304 with Acrylic Hood
- Temperature range ambient +5°C to 70°C
- Tray will be S.S. 304

Model No.	Size
OT-INS-61-MS/SS-9	9 Flask of 250 ml.
OT-INS-61-MS/SS-12	12 Flask of 250 ml.
OT-INS-61-MS/SS-16	16 Flask of 250 ml.

### Optional:

1. Digital RPM Indicator

MOC : MS-Mild steel duly powder coated SS-stainless steel kindly select the one required in model number.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@optics technology.in](mailto:info@optics technology.in),

[optics technology@gmail.com](mailto:optics technology@gmail.com),

Website : [www.optics technology.in](http://www.optics technology.in)