



## MUFFLE FURNACE CE



OT-MF-84 MS

M.S. BODY



OT-MF-84 SS

S.S. 304 BODY

Light Weight High Temp. Ceramic Fiber Wool Insulation. Outer Body Double Walled with M.S. Duly Powder Coated / Stainless Steel Body (S.S. 304). Heating Elements are of KHANTHAL "A-1" Wire, Temperature Control Unit Consist of Microprocessor based PID controller (Digital LED Display) fitted in front of Furnace, two pilot lamps & mains On/Off switch. Power plug to work on 220 Volts AC.

Maximum Temp. : 1000°C / 1200°C

Max. Working Temp. upto : 900°C / 1100°C

| Model No.        | Working Sizes<br>L x H x D         |
|------------------|------------------------------------|
| OT-MF-84 MS/SS/A | 4" x 4" x 9" (100 x 100 x 225 mm)  |
| OT-MF-84 MS/SS/B | 5" x 5" x 10" (125 x 125 x 250 mm) |
| OT-MF-84 MS/SS/C | 6" x 6" x 12" (150 x 150 x 300 mm) |
| OT-MF-84 MS/SS/D | 8" x 8" x 12" (200 x 200 x 300 mm) |
| OT-MF-84 MS/SS/E | 9" x 9" x 12" (225 x 225 x 300 mm) |
| OT-MF-84 MS/SS/F | 9" x 9" x 18" (225 x 225 x 450 mm) |

### Optional :

1. Programmable Controller with Ramp & Soak facility (Profile Controller)
2. Microprocessor based Touch Screen Display System
3. 21 CFR Software (For Non Editable Data)