



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

TEST TUBE WARMER CE

SPECIALLY FOR I.V.F. LABS



OT-TTW-229

TEST TUBE WARMER

Digital units include built-in temperature sensing probe Compact design saves valuable bench space.

Stainless Steel body construction ensures durability.

Digital units utilize a Digital Temperature Controller and easy-to-read LED display for accurate and reproducible temperature control.

Can hold 15ml test tube.

- 15 No. Tubes
- Imported sensor
- Temperature Range: ambient +5°C to 50°C.
- Temp. remains stable and locks at 37°C
- Heating blocks
- Accuracy $\pm 0.5^{\circ}\text{C}$

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in