



# OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED  
An ISO 10002:2018 CERTIFIED  
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010  
BS EN 12469:2000  
ISO 29463-5:2022

## WATER BATH HOLES (THERMOSTATICALLY CONTROLLED)

CE

Calibration is also Available



OT-WB-36-H/6-MS

THERMOSTATICALLY  
CONTROL (6 HOLES)

(Double walled) Inner chamber made of Stainless Steel 304 & outer made of S.S. 304 / Mild steel duly Powder Coated. The top has 75 mm holes with concentric rings. Temperature ranges upto 99°C. Accuracy  $\pm 5^\circ\text{C}$  in thermostatically Controlled fitted with Neon Lamp indicator, main switch and workable on 220 Volt, 50 Hz. A.C.

Model No.	Size	Holes
OT-WB-36 H/6-MS/SS	325x250x100mm	6 Hole of 75mm Dia
OT-WB-36 H/12-MS/SS	425x325x100mm	12 Hole of 75mm Dia
OT-WB-36 H/15-MS/SS	500x400x100mm	15 Hole of 75mm Dia

### Optional:

1. Microprocessor based PID temperature controller with an accuracy of  $\pm 0.5^\circ\text{C}$ .
2. Water level cutoff device.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in),

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com),

Website : [www.opticstechnology.in](http://www.opticstechnology.in)