



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

FLOCULATOR - JAR TESTER OR MULTIPLE SPINDLE STIRRER

CE



OT-FT-146

Fitted with D.C. MOTOR But to work on 220 Volts A.C. with Variable Speed from 20 R.P.M. to 100 R.P.M. approx. Jar Test Apparatus allows efficient and economical flocculation matter and collides in general and also has been specially designed for use in water treatments to correctly estimate the dosing of alum and such other coagulants for treatment of Water sewage and flocculation formation test in Water Treatment plants.

STIRRING CAPACITY

- (A) 2 Stirrers
- (B) 4 Stirrers
- (C) 6 Stirrers
- (D) 8 Stirrers

Available Capacity :

OT-FT-146/1 1 Ltr.

OT-FT-146/2 2 Ltr. etc.

Optional:

Mechanical Timer **OR** Digital Timer & Digital R.P.M. Indicator

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in