



SINCE 1978

# OPTICS TECHNOLOGY

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

**GLASS DISTILLATION (SINGLE) CE**

**WITH AUTOMATIC WATER LEVEL CUT OFF**



**OT-GD-S-121**

**SINGLE DISTILLATION**

The distillation apparatus consists of boiler with heating elements embedded in glass and fused in spiral type coil internally of the bottom and tapered round glass, suitable to work on 220 Volts 50 Cycles AC complete with metallic stand rings and clamps to hold boiler condenser electric board with switch and connecting cords. It is capable of producing pure distilled water. **Also the whole unit can be hanged on wall.**

**CAPACITY :** 4Ltr./hr. - Borosilicate Glass

**2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)**

**☎ 9811238150, 9811238151, 9811236960**

**E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)**

**[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)**

**Website : [www.opticstechnology.in](http://www.opticstechnology.in)**