



OPTICS TECHNOLOGY



SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

HEATING BLOCK / DRY BATH CE



OT-HB-227/1



OT-HB-227/2

(any Non Standard size)

This is a solid block heated unit with provision for different samples at a time. All the sample are heated at 150°C within $\pm 1\%$ accuracy with utmost uniformity of temperature.

These are fitted with microprocessor based PID Temperature controller

Heating block type 5°C

Temperature range from Ambient +5°C to 70°C $\pm 1^\circ\text{C}$

Operate on 230 volts, 50 Hz, single phase

Model No.	Available sizes for
OT-HB-227-A	6 samples
OT-HB-227-B	15 samples
OT-HB-227-C	20 samples
OT-HB-227-D	24 Samples at a time.
OT-HB-227-E	48 Samples at a time.
OT-HB-227-F	64 or any other non standard

AS PER REQUIREMENT WITH TUBE SAMPLES

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in