



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

KARL FISCHER TITRATOR (AUTOMATIC) (Microprocessor Based) CE



OT-MKFA-181

**MICROPROCESSOR
BASED KARL
FISCHER TITRATOR**

Microprocessor based technology, Motorized Dispensing & External Printer Attachment Facility.

Specifications:

Measuring Range	: 100 ppm to 100% of moisture
Moisture Content	: % age : (0.00 to 100%)
Calculation	: ppm : (100 to 999999)
Automatic	: mg H ₂ O : (0.00 to 999.99)
End Point Alarm	: Visual LED and Audible Alarm
Wait Time For End Point Alarm	: Adjustable (1 to 59 second)
Keyboard	: 24 keys soft touch, membrane type
Sample Storage	: Upto 20 latest samples
Display	: 8 digit bright red, seven segment LED
Clock	: In-built clock with date and time display
Stirrer	: In-built magnetic stirrer with adjustable speed

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in