



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

VACUUM OVEN (RECTANGULAR) CE



OT-RVO-137 MS



OT-RVO-137 SS

Inside chamber made of S.S. 304 Outside Mild Steel / S.S. 304 with Toughened Glass Door. It is provided with 1 or 2 No. of S.S. Trays.

Temp. Range Ambient Temp. +5°C to 130°C. Temp. Controlled by microprocessor based PID controller. Vacuum Pressure indication by vacuum gauge.

Model No.	W x D x H	Ltrs. (App.)
OT-RVO-137-MS/SS-A	200 x 200 x 200 mm (8" x 8" x 8")	8
OT-RVO-137-MS/SS-B	250 x 250 x 250 mm (10"x10"x10")	16
OT-RVO-137-MS/SS-C	300 x 300 x 300 mm (12"x12"x12")	28
OT-RVO-137-MS/SS-D	450 x 450 x 450 mm (18"x18"x18")	96
OT-RVO-137-MS/SS-E	600 x 450 x 450 mm (24"x18"x18")	128
OT-RVO-137-MS/SS-F	600 x 600 x 600 mm (24"x24"x24")	227

Optional:

1. Digital Vacuum Gauge.
2. Microprocessor Based Touch Screen Display System

Note: Without Vacuum Pump.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in