



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

BOX COMPRESSION TESTER

Microprocessor Based Model

CE



OT-BCT-192



OT-BCT-192

Capacity - 1000 Kgf

Test Speed - $\frac{1}{2}$ "/Min. and Location Speed 7.5"/Min.

Suitable to test for Boxes upto 100 x 100 x 100 cm.

Models Available:

- (a) Microprocessor based Digital Semi Automatic
- (b) Microprocessor based Digital Semi Automatic with Printer Port Facility
- (c) Microprocessor based Digital Fully Automatic with PC Control and operation
- (d) Microprocessor based Digital Fully Automatic with Touch Screen Monitor Control

Note: Printer will be arranged by the client

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in