



SINCE 1978

# OPTICS TECHNOLOGY



**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## COBB SIZING TESTER (UNIVERSAL MODEL) CE



The test is used for determining liquid absorptiveness or resilience of treated and untreated paper, board, fabrics and other sheet materials. The test consists of a test apparatus, its 10kg roller and blotting paper. The roller method is recommended to insure better consistency and reliability that can be obtained by manual blotting or rubbing the specimen.

### Specifications:

- Standard test area : 100 cm<sup>2</sup>
- Metal roller : 10 kg ± 0.5 kg

### Dimensions:

- Depth : 150 mm
- Height : 275 mm
- Width : 330 mm

### Features:

- One Universal model to test paper, Board and water proof board.
- Three graduations of height 10, 25 and 100 mm as specified in TAPPI 441 for Paper, Corrugated Board and Waterproof Board respectively.
- Adjustable clamp pressure mechanism for necessary required clamping pressure to prevent water leakage during clamping.
- Aluminium Block and steel Roller for rust proof and longer life.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)