



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2013 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DROP TESTER

CE



OT-DT-197



OT-DT-197

Drop tester manual / electrically operated auto model for determination of transport worthiness of packed containers.

Features:

- Flap Opening Type and Manual height adjustment and locking system
- It is highly consistent equipment for findings repeatable and accurate results.
- The testing instruments supplied with conformance certificate.
- Drop Height adjustable clamp
- Strong base plate with rugged structure

Specification:

Max. Drop height - 1800 mm & minimum drop height 700 mm.

Interval = 1 mm

Maximum Box Weight at: 50 Kgs.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in