



# OPTICS TECHNOLOGY

SINCE 1978

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED  
An ISO 10002:2018 CERTIFIED  
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010  
BS EN 12469:2000  
ISO 29463-5:2022

## PUNCTURE RESISTANT TESTER

CE



OT-PRT-193

The puncture resistance Tester is used for determine the puncture of corrugated board by measuring of energy required to puncture it. Interchangeable weight supplied to cover four energy ranges i.e. in new meter 0-5, 0-10, 0-20 & 0-36.

### Applications for:

Paper Industries and Packaging Industries

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)