



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED
An ISO 10002:2018 CERTIFIED
An ISO 13485:2016 CERTIFIED

IEC 61010-1:2010
BS EN 12469:2000
ISO 29463-5:2022

AUTOMATIC WHITENESS METER

CE

(Table Top Model) Mains Operated



OT-AWM-27 TT

TABLE TOP

Technical Details:

Measurement range	: 0 ~ 100%
Whiteness formula	: blue light whiteness
Light source	: LED light source
Display	: LCD display
Measuring the aperture	: 20
Minimum	: 0.1%
Measurement repeatability	: 0.2%
Power supply	: DC 1.5V x 5 AA alkaline dry

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in