



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

BURSTING STRENGTH TESTER

(ANALOGUE MODEL)



OT-BST-194/1

The equipment is used for accurate measurement of Bursting Strength of Paper, Paper Boards, Corrugated Boards, Solid Fibre Board, Filter Cloth/ Industrial Fabric etc. and is manufactured as per national and various International standards, viz., TAPPI, ASTM, FEFCO etc.

- Fully enclosed model.
- Bursting Strength displayed on gauges for greater accuracy.
- Accurate, reliable and trouble free equipment.
- Available two pressure gauges.
- Change over from one gauge to the other by a simple switch.
- Machine available in various specifications to test Paper, Paper Board, Corrugated Board, Solid Fibre board, Filter Cloth / Industrial Fabric etc.
- Range 0-7 Kg in one gauge and 0-35 Kg in another

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in