



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

POCKET CONDUCTIVITY METER



Measuring range	0 – 19990 μ S/cm	Temp. : 0 - 50°C
Resolution	10 μ S/cm	Temp. : 0.1°C
Accuracy	+/-2% F.S.	Temp. : +/-1°C

Single point manual calibration

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in