



DEMINERALISER / DE-IONIZER

CE



OT-DM-126

It consist of an inbuilt prefilter & very high exchanging capacity non-corrosive CATION & ANION Resin, columns prefilled with imported Resins for efficient results having treated water parameters. Fitted with a sturdy PVC multi control valve & can be connected direct to water taps. Fitted with analog / digital conductivity meter. Supplied with one plastic chemical proof regeneration tank & the complete unit is housed in a painted M.S. trolley for easy mobility.

FEATURES:

- Digital Display
- Non corrosive Regeneration Tank
- In Built Pre filter
- High capacity Cation and Anion exchange resins.
- Upflow treatment.
- Easily portable - metal stand and trolley
- Can be connected direct to Water Tap.
- Simple Periodic regeneration - negligible operating costs.
- Digital model is fitted with digital conductivity meter.
- Treated water will have the conductivity less than 10 ms/cm and pH of 7.5 - 9

• Available in **Double Bed**, **Four Bed** and **Mixed bed** models

Capacity:

- 50 Ltrs. / Hr.
- 70 Ltrs. / Hr.
- 100 Ltrs. / Hr.
- 200 Ltrs. / Hr.
- 500 Ltrs. / Hr.