



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DEMINERALISER / DE-IONIZER

CE



OT-DM-126

It consist of an inbuilt prefilter & very high exchanging capacity non-corrosive CATION & ANION Resin, columns prefilled with imported Resins for efficient results having treated water parameters. Fitted with a sturdy PVC multi control valve & can be connected direct to water taps. Fitted with analog / digital conductivity meter. Supplied with one plastic chemical proof regeneration tank & the complete unit is housed in a painted M.S. trolley for easy mobility.

FEATURES:

- Digital Display
- Non corrosive Regeneration Tank
- In Built Pre filter
- High capacity Cation and Anion exchange resins.
- Upflow treatment.
- Easily portable - metal stand and trolley
- Can be connected direct to Water Tap.
- Simple Periodic regeneration - negligible operating costs.
- Digital model is fitted with digital conductivity meter.
- Treated water will have the conductivity less than 10 ms/cm and pH of 7.5 - 9

• Available in **Double Bed**, **Four Bed** and **Mixed bed** models

Capacity : 50 70 100 200 500 Ltrs. / Hr.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)
☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in