



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL CONDUCTIVITY AND TDS METER COMBINED

(Table Top) CE



OT-DCTDS-109



OT-DCTDS-109-M

DETAILS	CONDUCTIVITY	TDS
Display	3½ Digit bright red LED display	
Resolution	0.1µS/cm in lowest range	0.1ppm in lowest range
Accuracy	± 0.5% FS ± 1 digit	
Ranges	0 to 200µS/cm	0 to 200ppm
	0 to 2mS/cm	0 to 2ppt
	0 to 20mS/cm	0 to 20ppt
	0 to 200mS/cm	0 to 200ppt
	0 to 1000mS/cm	0 to 1000ppt

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in