



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL CONDUCTIVITY METER (Table Top) CE



OT-COM-106-D

Ranges

0 to 0.2 mS/cm	(0 to 200 μ S/cm)
0 to 2 mS/cm	(0 to 2000 μ S/cm)
0 to 20 mS/cm	(0 to 20000 μ S/cm)
0 to 200 mS/cm	(0 to 200000 μ S/cm)
0 to 1000 mS/cm	(0 to 1000000 μ S/cm)

Accuracy

$\pm 1\%$ FS ± 1 Digit

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)
☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in