



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DIGITAL POTENTIOMETER

(Table Top)

CE



OT-DP-110

Features:

- » Compact Designed
- » 3½ Digit LED Display
- » Highly Stable & Accurate
- » Rugged & Economical
- » 4 Electrodes for Different Titrations

Specifications:

Range	: 0 – ± 199.9 mV, 0 – ± 1999 mV
Resolution	: 0.1 mV
Repeatability	: ± 1 mV
Accuracy	: ± 1 mV, ± 1 digit
Power	: 230 V ± 10% AC, 50 Hz.
Temperature	: 10°C - 45°C
Input Impedance	: >10 ¹² ohms

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in