



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DISTILLATION APPARATUS

CE

TABLE TOP MODEL (HEAVY DUTY)



OT-DATP-123



Water level cutoff



ALSO USED FOR
HERBAL EXTRACT

OT-DATP-123

- Double Walled Construction.
- Inner & Outer Made of S.S. 304.
- Deluxe Model Heavy Duty with Argon Welding.
- The gap in between being filled with glass wool.
- The condenser & TOP are also made of S.S.
- To work on 220/230 volts.

Model No.	Capacity
OT-DATP-123/A	4 Ltr./hr
OT-DATP-123/B	6 Ltr./hr
OT-DATP-123/C	8 Ltr./hr
OT-DATP-123/D	10 Ltr./hr
OT-DATP-123/E	20 Ltr./hr
OT-DATP-123/F	25 Ltr./hr
OT-DATP-123/G	40 Ltr./hr
OT-DATP-123/H	50 Ltr./hr
OT-DATP-123/I	80 Ltr./hr
OT-DATP-123/J	100 Ltr./hr

Optional:

1. Water level cutoff device (For Single Phase)
2. Water level cutoff device (For Three Phase)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in