



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

ELECTRO MAGNETIC SIEVE SHAKER (NOISELESS) CE



OT-EMSS-115

The instrument is powered by an electromagnetic drive which has no rotating parts, making it maintenance free and extremely quiet in operation.

The unit has been designed for user friendly operation. The top plate of the Sieve shaker has clearly marked calibrated lines to aid visual inspection of the amplitude.

- Sieving with 3-D effect
- For sieve 200mm (8")
- Measuring range 20mm to 10mm
- Digital setting and control of sieving parameters
- Reproducible and globally comparable sieving results
- Low noise, maintenance - free

Optional: • Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in